


INSTRUCTIONS FOR USE

Dupixent 300 mg Single-use pre-filled pen

The Instructions for Use contain information explaining how to inject Dupixent.

Read the instructions before using the pen. **Do not inject yourself or anyone else until you have been shown how to inject Dupixent.** Your medical team can show you or your caregiver how to prepare and inject a dose of Dupixent before you try to do it yourself for the first time. Keep these instructions for future use. Call your medical team if you have questions. **The Dupixent pre-filled pen is intended for use in adults and children from the age of 6 years and above.** The Dupixent pre-filled pen is a single-use device. It contains 300 mg of Dupixent for injection under the skin (subcutaneous injection).

The parts of the Dupixent pre-filled pen are displayed below:



Important information

- Read all of the instructions carefully before using Dupixent pre-filled pen.
- Ask your healthcare provider how often you will need to inject the medicine.
- **In children aged 12 years and older**, it is recommended that Dupixent be administered by or under the supervision of an adult.
- **In children 6 years to 12 years of age**, the Dupixent injection must be administered by a caregiver.
- Choose a different injection site for each injection.
- **Do not** press or touch the yellow needle cover with your fingers.
- **Do not** inject through clothes.
- **Do not** remove the green cap until just before the injection.
- **Do not** try to put the green cap back on the pen.
- Throw away the used pen immediately after use.
- **Do not** re-use a Dupixent pre-filled pen.

Storing Dupixent

- Before use, store pens in the refrigerator at a temperature between 2°C-8°C.
- Store pens in the original package in order to protect from light.
- If necessary, Dupixent pens can be stored at room temperature (up to 25°C) for **up to 14 days**.
- **Do not** store Dupixent pens at a room temperature above 25°C.
- After removing a Dupixent pen from the refrigerator, use it within 14 days or throw it away.
- **Never** shake the pen.
- **Do not** heat the pen.
- **Do not** freeze the pen.
- **Do not** put the pen into direct sunlight.
- **Keep the Dupixent pens and other medicines out of the reach and sight of children.**

9094017-R

A. GET READY TO INJECT

A1: Prepare supplies

Choose a clean and flat work surface. Ensure you have the following items:
the Dupixent pen
*an alcohol wipe
*a cotton ball or gauze pad
*a puncture-resistant container (see Step D)
**These items are not included in the product package*

A2: Check the pen

Do not use the Dupixent pen if it has been damaged. **Do not** use the Dupixent pen if the green cap is missing or not securely attached. **Do not** use the Dupixent pen if the window is yellow before use.

A3: Look at the label

Check the preparation and the dose

- Make sure you have the correct preparation and the correct dose.



- Check the expiration date.
- ⚠ **Do not** use the Dupixent pen if the expiration date has passed.



A4: Check the medicine



- Look at the medicine through the window. It should be clear and colorless to pale yellow.
- Note: You may see an air bubble, this is normal.
- ⚠ **Do not** use the pen if the liquid is discolored or cloudy, or if it contains visible flakes or particles.


A5: Wait 45 minutes



- Place the Dupixent pen on a flat surface and let it warm up naturally at room temperature, below 25°C, for at least 45 minutes.
- ⚠ **Do not** heat the Dupixent pen in a microwave, hot water or direct sunlight.

B. CHOOSE AND PREPARE YOUR INJECTION SITE

B1: Wash your hands well with soap and water




B2: Choose an injection site

- **Thigh**
- **Abdomen**, except for 5 cm around your navel.
- A caregiver can also inject in the outer area of your **upper arm**.
- Choose a different injection site each time you inject Dupixent.
- ⚠ **Do not** inject into skin that is tender, damaged, has bruises or scars, or into areas with visible veins.
- ⚠ **Do not** inject through clothes.


B3: Prepare the injection site

- Clean the injection site with an alcohol wipe.
- Let the skin dry before injecting.
- ⚠ **Do not** touch the injection site again or blow on it before the injection.



C. GIVE THE INJECTION

C1: Remove green cap




- Remove the green cap by pulling it straight off, as shown in the illustration. **Do not** twist the green cap off.
- **Do not** remove the green cap until you are ready to inject.
- **Do not** touch or press the yellow needle cover with your fingers. The needle is inside.
- ⚠ **Do not** put the green cap back on the pen after you have removed it.


C2: Pinch the skin and place

Pinch the skin before and during the injection.

- Hold the Dupixent pen so that you can see the window. Place the yellow needle cover on the skin.
- Place the yellow needle cover on the skin at approximately a 90-degree angle.
- Pinching of the skin is not needed for adults and children aged 12 years and older.



C3: Press down → Look at the window until it turns completely yellow → Then count to 5



Press and hold the Dupixent pen firmly against the skin **until you cannot see the yellow needle cover**.

- There will be a “click” when the injection starts, **and then**
- The window will start to turn yellow.
- **Keep pressing** the Dupixent pen against the skin.


Keep pressing the Dupixent pen against the skin and **watch the window**:

- The window will turn completely yellow, **and then**
- You will hear a **second “click”**.
- **Keep pressing** the Dupixent pen against the skin.

Keep pressing the Dupixent pen against the skin **and count to 5 to make sure you get the full dose**.

Pinching of the skin is not needed for adults and children aged 12 years and older.

C4: Remove



- After you have completed the injection, pull straight up to remove the Dupixent pen from the skin. The yellow needle cover will cover the needle.
- If you see blood at the injection site, gently wipe with a cotton ball or gauze pad.
- ⚠ **Do not** rub the skin after the injection.

If the window does not turn completely yellow, or if it looks like medicine is still coming out of the pen, you may not have received the full dose. Dispose of the pen and contact your medical team immediately.


⚠ **Do not** administer a second dose without speaking to your medical team.

D. DISPOSAL OF A USED DUPIXENT PEN

How to dispose of the Dupixent pen

Dispose of used Dupixent pens and green caps in a puncture-resistant container, immediately after use.

- ⚠ **Do not** dispose of used Dupixent pens and green caps in the household waste.



- ⚠ **Do not put the green cap back on the pen.** Keep the container out of the reach and sight of children.