

Anticoagulant Citrate Dextrose Solution, USP (ACD) Formula A

Anticoagulant Citrate Dextrose Solution is sterile, nonpyrogenic and for use with Cytapheresis Device only. It is not for direct intravenous infusion.

The plastic Container is fabricated from a specially formulated polyvinyl chloride. Water can permeate from inside the container into the overwrap in amounts insufficient to affect the solution significantly. Solution in contact with the plastic container can leach out certain of its chemical components in very small amounts within the expiration period, e.g., di-2-ethylhexyl phthalate (DEHP), up to 5 parts per million. However, the safety of the plastic has been confirmed in tests in animals, according to USP biological standards for plastic container as well as by tissue culture toxicity studies.

Do not use unless solution is clear and no leaks are detected.

Single use container; discard unused portion.

Do not remove unit from overwrap until ready for use. Use all units promptly when pouch is opened.

The overwrap is a moisture barrier. The inner bag maintains the sterility of the product.

After removing overwrap, check for minute leaks by squeezing inner bag firmly. If leaks are found, discard solution as sterility may be impaired.

Store at room temperature. Avoid storage at excessive heat.
Protect from freezing.

Drug Reg. No.: 111 96 27598 00

Teva Medical Ltd,
P.O.Box 2, Ashdod